


PulsarTM

Self-Expanding Stent System

Thin struts, low chronic outward force (COF).

36-month results confirm the excellent
12- and 24-month outcomes of the PulsarTM stent.




Thinner struts and low COF make a difference

140 µm thin struts – thinner than the leading brands¹


Chronic outward force (COF)

Chronic outward force (COF) is the force exerted on the vessel wall by a self-expanding (SE) stent to achieve its preset diameter.



SE Stent 1 mm oversized in vessel

Minimal metal burden¹




*Bare metal platform for Zilver PTX drug eluting stent

Thinner struts for low chronic outward force (COF)²

With low COF the natural response of the body is much lower when compared to high COF stents^{**}:

- less inflammation³
- less neointimal hyperplasia³
- lower risk of restenosis³


^{**}As demonstrated in pre-clinical studies



Sufficient radial force for long-term vessel support, even in calcified lesions

The following figure shows a fluoroscopic image after stent implantation (6 × 80 mm Pulsar™-18 in 2011) in a severely calcified lesion, with lack of full expansion of the stent in an area of focal calcification, and an angiographic residual stenosis of less than 30 % (circle; with a high COF stent the stress exerted on the vessel wall would be significant, potentially leading to neo-intimal hyperplasia).

A fluoroscopic image that was obtained 5 years later when the patient came back for treatment of an ipsilateral iliac lesion shows full expansion of the stent, and angiographically a full patency of the stented segment. This figure illustrates nicely the chronic impact of the Pulsar™ stent and its COF on long-term clinical outcome. With a constant low chronic outward force applied to the vessel, patency can be achieved and maintained over a long term follow-up.



After the treatment 2011

2016

(Courtesy of Prof. van den Berg)

Stent strut thickness in perspective¹




Benefits of low COF proven in animal studies

Stents with higher COF show more neointimal hyperplasia^{**}


Pulsar™ stent causes less neointimal hyperplasia^{4,**}

Pre-clinical study results indicate that the low COF Pulsar™ stent causes less neointimal hyperplasia at 90 days compared to the high COF Lifestent stent.



Neointimal area vs. COF^{3,**}

Stents with higher COF show more neointimal area.



Thinner struts and lower COF make a difference:^{**}


- Lower risk of restenosis³
- Reduced vessel injury and inflammation³
- Faster endothelialization^{5,6}

^{**}As demonstrated in pre-clinical studies

Pre-clinical study results showed significantly smaller area stenosis, lower injury and inflammation scores for the Pulsar™ stent compared to the high COF stent Lifestent.^{**}

Benefits of low COF shown in clinical study

BIOFLEX-I⁷: COF correlation to peak systolic velocity ratio (PSVR)




- At 1 month COF shows minimal impact on blood flow
- At 12 months blood flow worsens with the increasing COF


Conclusion
 In the long term, vessel reaction on high COF might lead to a higher risk of restenosis.

Risk factors for increased PSVR		
MODEL		SIG.
Radial outward force		0.024

COF was one of the most significant predictors for deterioration of PSVR.



Bending
 Peak-to-valley design and S-articulating connecting bars provide multi-directional flexibility and avoid fish scaling in mobile vessel architecture.¹




Low COF – Less is more.



Clinical evidence


Low COF, thin struts – proven in extensive clinical program

Consistent 12-month outcomes across different lesion length			
STUDY	PRODUCT	ALL ¹	PP ¹¹
BIOFLEX PEACE ⁸	Pulsar™ (stent only)	8.2 cm	84.7 %
4EVER ⁹	Pulsar™	7.1 cm	81.4 %
BIOFLEX-I ¹⁰	Pulsar™	8.2 cm	70.6 %
TASC D II ¹¹	Pulsar™	18.2 cm	85.4 %
TASC D ¹²	Pulsar™	24.5 cm	77.0 %
PEACE ¹³	Pulsar™	11.2 cm	79.5 %




Excellent outcomes after 2 years

24-month data in perspective			
STUDY	PRODUCT	ALL ¹	PP ¹¹
BIOFLEX PEACE ⁸	Pulsar™ (stent only)	8.2 cm	78.4 %
SUPERB ¹⁴	Supera	7.8 cm	n/a
4EVER ¹⁵	Pulsar™	7.1 cm	72.3 %
BIOFLEX-I ¹⁰	Pulsar™	8.2 cm	n/a
STROLL ¹⁶	S.M.A.R.T. Control	7.7 cm	74.9 %
RESILIENT ¹⁷	Lifestent	7.0 cm	n/a
ZILVER PTX ¹⁸	Zilver BMS provisional	6.3 cm	65.8 %
DURABILITY II ¹⁹	EverFlex	10.9 cm	66.1 %



¹Average Lesion Length ¹¹Primary Patency ¹¹¹Freedom from Target Lesion Revascularization
 Results from different trials are not directly comparable. Differences in outcomes may be the result of differences in protocol design, patient populations or other factors. Astron Pulsar™, Pulsar™-18, Pulsar™-18 T3 and Pulsar™-35 have equivalent stent platforms, therefore the clinical results are valid for the Pulsar™ range.

3-year outcomes in perspective



After 3-years fTLR is still in line with competitor bare metal stent studies.

Pulsar™ stent’s thin strut and low COF efficiency are proven throughout stable outcomes over the long term.

Key studies

Outcomes of key Teleflex studies

	TASC D ¹²	4EVER ⁹	BIOFLEX PEACE ⁸	BIOFLEX-I ¹⁰
STUDY	LONG AND OCCLUDED LESIONS IN HIGH RISK POPULATION	ENDOASCULAR TREATMENT WITH 4F DEVICES	STENT ONLY GROUP	IDE-TRIAL
Selected Demographics	– 100 % CLI ³ patients – 72.7 % Diabetics – 100 % Occlusion	– 16.7 % CLI ³ patients – 35.8 % Diabetics – 20.8 % Occlusions	– 15.8 % CLI ³ patients – 40.0 % Diabetics – 76.7 % Smokers	– 40.7 % Diabetics – 30.1 % Total occlusion
Patient Number	22	120	60/160#	302
Average Lesion Length	24.5 cm	7.1 cm	8.2 cm	8.2 cm
	TASC D ¹²	4EVER ⁹	BIOFLEX PEACE ⁸	BIOFLEX-I ¹⁰
12-MONTH RESULTS				
Primary Patency	77.0 %	81.4 %	84.7 %	70.6 %
FTLR ¹¹	86.0 %	89.3 %	89.3 %	87.6 %
Additional Outcomes	– Limb Salvage: 100 %	– Access site complication: 0.9 % (excludes patients on coumarin) – Mean Compression Time: 8.1 minutes	– Oversizing of 0.8 mm – 81.7 % Improvement of Rutherford classification	– Study results are similar to the results of other US IDE trials
	TASC D ¹²	4EVER ¹⁵	BIOFLEX PEACE ⁸	BIOFLEX-I ¹⁰
24-MONTH RESULTS				
Primary Patency	n/a	72.3 %	78.4 %	n/a
FTLR ¹¹	n/a	82.7 %	89.3 %	81.0 %
36-MONTH RESULTS				
FTLR ¹¹	n/a	n/a	n/a	78.2 %
Key Messages	– Pulsar™ stent has shown excellent outcomes for long occluded lesions (100 % TASC D lesions) in TASC D study. – Pulsar™ stent has shown to reduce the risk of foot loss – Pulsar™ stent has shown its efficacy in a high risk patients population (Mean Rutherford class by 100 % CLI patients)	– Improved acute outcomes of 4F vs. 6F ¹ – Less access site complication rates may permit ambulatory treatment – Shorter compression time with 4F compared to 6F – 30.8 % Calcification: No significant difference between calcified vs. non-calcified lesions at both 12 m and 24 m	– Pulsar™ stent is very efficient under All-comers condition – Pulsar™ stent shows excellent results in mixed population	– FDA approved – Pulsar™ stent systems are the only US SFA stents that are compatible with 4F sheaths ¹¹ – The study results are similar to the results of other competitor stents for use in patients with femoropopliteal lesion – 36 m results confirm Pulsar™ stent safety and efficiency on the long term

³CLI–Critical Limb Ischemia ¹Less access site complications ¹¹when compared to leading US brands ¹¹¹Freedom from Target Lesion Revascularization

References

- 1 Data on file. 6.0 mm diameters.
- 2 Data on file. 6.0 mm diameters. Supera stent not possible to test due to its design and applied test method.
- 3 Zhao HQ Late stent expansion and neointimal proliferation of oversized nitinol stents in peripheral arteries. Cardiovasc. Interv. Radiol. 2009; 32(4): 720-6.
- 4 Funovics M. Correlation between chronic outward force (COF) and neointimal hyperplasia in self-expanding nitinol stents in swine in clinically relevant oversizing ranges. Presented at: LINC, Jan 26, 2017; Leipzig, Germany.
- 5 Koskinas C. Role of endothelial shear stress in stent restenosis and thrombosis: pathophysiologic mechanisms and implications for clinical translation. JACC 2012 10;59(15):1337-49.
- 6 Koppa T. Thrombogenicity and early vascular healing response in metallic biodegradable polymer-based and fully bioabsorbable drug-eluting stents. Circ Cardiovasc Interv. 2015 8(6):e002427.
- 7 Burket M BIOFLEX-I Increased outward force of SE-BMS in the SFA. Presented at: CX, 16 Apr, 2019. London, UK.
- 8 Lichtenberg et al. Effectiveness of the Pulsar-18 self-expanding stent with optional drug-coated balloon angioplasty in the treatment of femoropopliteal lesions–the BIOFLEX PEACE All-Comers Registry.Vasa (2019), 1-9. doi_10.10240301- 1526a000785.
- 9 Bosiers M et al. 4-French –compatible endovascular material is safe & effective in the treatment of femoropopliteal occlusive disease: Results of the 4EVER Trial. ENDOVASC THER 2013; 20: 746-756.
- 10 BIOFLEX-I Pulsar 2018 Post Approval Clinical Report_Final_36m.
- 11 Lichtenberg M et al. Evaluation of the 4-French Pulsar-18 Self-expanding Nitinol Stent in Long Femoropopliteal Lesions. Clinical Medicine Insights: Cardiology. 2014; 8(S2): 37 -42; 12. Lichtenberg M. Superficial Femoral Artery TASC D registry: 12-month effectiveness analysis of the Pulsar-18 SE nitinol stent in patients with critical limb ischemia. J Cardiovasc Surg (Torino). 2013 ; 54(4):433-9.
- 13 Lichtenberg M. BIOFLEX PEACE registry: 12 and 24 month results. Presented at:LINC 2018. Leipzig, Germany.
- 14 Supera IFU, EL2100430 (2016-03-23).
- 15 Bosiers M. 4EVER 24 month results: long-term results of 4F Pulsar stent in femoropopliteal lesions. Presented at: CIRSE 2013;Barcelona, Spain.
- 16 Bunte M et al. in STROLL Catheterization and Cardiovascular Interventions 2018; 92:106-114.
- 17 Laird J et al. RESILIENT SFA nitinol stenting. JET 2012;19:1-9.
- 18 Dake M et al. Durable clinical effectiveness with paclitaxel-eluting stents in the femoropopliteal artery: 5-year results of the Zilver PTX randomized trial. Am Heart Assoc 133.15 (2016); 1472-1483. doi:10.1161/CIRCULATIONAHA.115.016900.
- 19 Rocha-Singh et al. DURABILITY II Three-Year Follow-up. Catheterization and Cardiovascular Interventions 2015; 86:164-170.

Refer to IFU

Refer to QR code for a copy of the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.



About Teleflex

As a global provider of medical technologies, Teleflex is driven by our purpose to improve the health and quality of people's lives. Through our vision to become the most trusted partner in healthcare, we offer a diverse portfolio with solutions in the therapy areas of anesthesia, emergency medicine, interventional cardiology and radiology, surgical, vascular access, and urology. We believe that the potential of great people, purpose driven innovation, and world-class products can shape the future direction of healthcare.

Teleflex is the home of Arrow™, Barrigel™, Deknatel™, LMA™, Pilling™, QuikClot™, Rüschi™, UroLift™ and Weck™ – trusted brands united by a common sense of purpose.

At Teleflex, we are empowering the future of healthcare. For more information, please visit teleflex.com.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Teleflex, the Teleflex logo, and Pulsar are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries. All other names are the trademarks or registered trademarks of their respective owners.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative.

© 2025 Teleflex Incorporated. All rights reserved.

10 25 PDF