





Improved
crossability



Excellent
deliverability



Low profile, wide
range of sizes



Technical data/
ordering info

Teleflex[™]
Empowering the future of healthcare

Passeo[™] -35 Xeo[™]

PTA Balloon Catheter

Low profile PTA balloon with excellent deliverability
in a wide range of sizes.






Passeo™-35 Xeo™ Balloon

Low profile PTA balloon with excellent deliverability in a wide range of sizes.


Improved crossability

Up to 0.16 mm lower tip entry profile compared to leading competitors.¹




Advanced tip profile

The low tip entry profile and smooth tapered tip of Passeo™-35 Xeo™ Balloon facilitates crossing.



Lowest crossing profile


Up to 0.22 mm lower crossing profile compared to leading competitors.¹






Excellent deliverability

The Passeo™-35 Xeo™ Balloon Catheter offers excellent trackability and pushability compared to competitive 0.035" PTA balloon catheters.²



Low profile, wide range of sizes


Passeo™-35 Xeo™ Balloon features a low profile design with a full 5F compatibility up to 7.0 mm balloon diameter and 250 mm balloon length, 6F up to 10.0 mm and 7F up to 12.0 mm diameter.





Passeo™ -35 Xeo™ Balloon

Low profile PTA balloon designed to cross



Optimal deflation times
due to size specific shaft design

5F compatibility
Up to 7.0 × 250 mm

250 mm
Balloon length up to 250 mm

ø 12.0 mm
Balloon diameter up to 12.0 mm

Low entry profile and smooth taper
for better crossing

Low crossing profile
to ease lesion crossing

Hydrophobic coating
Patchwork design

Radial approach
Usable length up to 170 cm


New hub design
with improved information readability

5/120
Passeo-35 Xeo
1018426-03
9NP atm 19RBP atm
0.035"



Passeo™-35 Xeo™ Balloon

Vascular
Intervention
Peripheral



The Passeo™-35 Xeo™ Balloon is indicated to dilate stenosis in the iliac*, femoral, popliteal and infra-popliteal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. Passeo™-35 Xeo™ Balloon is also recommended for post-dilation of balloon expandable and self-expanding stents in the peripheral vasculature.**

Technical data

BALLOON CATHETER	
Catheter type	OTW
Recommended guide wire	0.035"
Tip	Low entry profile, colored
Balloon material	SCP (Semi-Crystalline Polymer), controlled compliance
Balloon coating	Hydrophobic patchwork coating
Balloon markers	2 swaged markers
Sizes	ø 3.0–12.0 mm; L: 20–250 mm
Shaft	5.1-5.4F, dual-lumen, hydrophobic coating
Usable length	90, 130 and 170 cm
Guide wire lumen	Hydrophobic coating

Ordering Information

BALLOON Ø		NP [‡]	RBP [‡]	CATHETER LENGTH 90 CM, BALLOON LENGTH				
				20 mm	40 mm	60 mm	80 mm	100 mm
5F	3.0 mm	9	21	428777	428786	428795	428804	-
	4.0 mm	9	20	428778	428787	428796	428805	428814
	5.0 mm	9	19	428779	428788	428797	428806	428815
	6.0 mm	9	16	428780	428789	428798	428807	428816
	7.0 mm	9	15	428781	428790	428799	428808	428817
6F	8.0 mm	9	14	428782	428791	428800	428809	428818
	9.0 mm	9	12	428783	428792	428801	428810	428819
	10.0 mm	8	11	428784	428793	428802	428811	428820
7F	12.0 mm	8	10	-	428794	428803	428812	428821

BALLOON Ø		NP [‡]	RBP [‡]	CATHETER LENGTH 90 CM, BALLOON LENGTH				
				120 mm	150 mm	170 mm	200 mm	250 mm
5F	3.0 mm	9	21	-	-	-	-	-
	4.0 mm	9	20	428823	-	-	-	-
	5.0 mm	9	19	428824	428833	-	428843	428848
	6.0 mm	9	16	428825	428834	-	428844	428849
	7.0 mm	9	15	-	-	-	-	-
6F	8.0 mm	9	14	-	-	-	-	-
	9.0 mm	9	12	-	-	-	-	-
	10.0 mm	8	11	-	-	-	-	-
7F	12.0 mm	8	10	-	-	-	-	-

BALLOON Ø		NP [‡]	RBP [‡]	CATHETER LENGTH 130 CM, BALLOON LENGTH				
				20 mm	40 mm	60 mm	80 mm	100 mm
5F	3.0 mm	9	21	428851	428860	428869	428878	428887
	4.0 mm	9	20	428852	428861	428870	428879	428888
	5.0 mm	9	19	428853	428862	428871	428880	428889
	6.0 mm	9	16	428854	428863	428872	428881	428890
	7.0 mm	9	15	428855	428864	428873	428882	428891
6F	8.0 mm	9	14	428856	428865	428874	428883	428892
	9.0 mm	9	12	428857	428866	428875	428884	428893
	10.0 mm	8	11	428858	428867	428876	428885	428894
7F	12.0 mm	8	10	-	428868	428877	428886	428895

BALLOON Ø		NP [‡]	RBP [‡]	CATHETER LENGTH 130 CM, BALLOON LENGTH				
				120 mm	150 mm	170 mm	200 mm	250 mm
5F	3.0 mm	9	21	428896	428905	428910	428915	428920
	4.0 mm	9	20	428897	428906	428911	428916	428921
	5.0 mm	9	19	428898	428907	428912	428917	428922
	6.0 mm	9	16	428899	428908	428913	428918	428923
	7.0 mm	9	15	428900	428909	428914	428919	428924
6F	8.0 mm	9	14	428901	-	-	-	-
	9.0 mm	9	12	428902	-	-	-	-
	10.0 mm	8	11	428903	-	-	-	-
7F	12.0 mm	8	10	428904	-	-	-	-

BALLOON Ø		NP [‡]	RBP [‡]	CATHETER LENGTH 170 CM, BALLOON LENGTH				
				20 mm	40 mm	60 mm	80 mm	100 mm
5F	4.0 mm	9	20	-	428935	428944	428953	-
	5.0 mm	9	19	428927	428936	428945	428954	-
	6.0 mm	9	16	428928	428937	428946	428955	-
	7.0 mm	9	15	428929	428938	428947	428956	-
6F	8.0 mm	9	14	428930	428939	428948	428957	-
	9.0 mm	9	12	-	428940	428949	428958	-
	10.0 mm	8	11	-	428941	428950	428959	-

BALLOON Ø		NP [‡]	RBP [‡]	CATHETER LENGTH 170 CM, BALLOON LENGTH				
				120 mm	150 mm	170 mm	200 mm	250 mm
5F	4.0 mm	9	20	428971	428980	-	428990	428995
	5.0 mm	9	19	428972	428981	-	428991	428996
	6.0 mm	9	16	428973	428982	-	428992	428997
	7.0 mm	9	15	-	-	-	-	-
6F	8.0 mm	9	14	-	-	-	-	-
	9.0 mm	9	12	-	-	-	-	-
	10.0 mm	8	11	-	-	-	-	-

[‡] NP=Nominal Pressure and RBP=Rated Burst Pressure in atm (1 atm=1.013 bar)
* Note for Australia: Passeo™-35 Xeo™ Balloon not approved by Therapeutic Goods Administration for use in the common iliac arteries.
** Indication as per IFU.

References:

1 Data on file. Compared to leading competitors 6.0 mm diameter balloon size. IIB(P) 39/2017.
2 Data on file. Compared to leading competitors 6.0 x 200 mm balloon size.

Teleflex, the Teleflex logo, Passeo and Xeo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries. All other trademarks marked with a ™ are the property of their respective owners and are solely used for identification purposes and do not imply any affiliation, endorsement, or ownership by Teleflex Incorporated or its affiliates. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative.

© 2026 Teleflex Incorporated. All rights reserved.

446265-EN · REV F · 01 26 PDF · DV

Distributed by:

Teleflex Headquarters International, Ireland · Teleflex Medical Europe Ltd. · IDA Business & Technology Park
Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com
United Kingdom Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com
South Africa Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com

teleflex.com

Teleflex™
Empowering the future of healthcare